

**NEW!**

# Rechargeable LED Head Light

## Operating LED Head Light (KS-05)



### Strong Point

- Brightness is more than 50,000 Lux
- Semi-permanent life with usage of LED lamp (100,000hours)
- No UV radiation
- Suitable for operating because of no heat from the lamp
- Use Li-ion battery that can be recharged more than 800 times
- 3-Hours of continual use with one recharge
- Delicate design without any exposed wires
- Excellent heatsink design and economical price
- The head band is made of soft natural leather
- Sliding type of both eye loupe equipable (option)
- Light Spot (9cm at 25cm)
- Stepless light intensity control (0~100%)
- Lens assembly : scientifically arranged three optical lens

### Model

Model	Light color	Type of battery	Accessories
KS-05	White	Battery integrated type	①×2, ②, ③
		Belt clip battery pack type	④, ⑥, ⑦
		Total type	①×2, ②, ③, ④, ⑥, ⑦



### Accessories



① Integrated Battery(2,200mAh)

② Battery charging holder

③ Charging Adapter(Integrated type)

④ Belt clip Battery pack(4400mAh)

⑤ Belt clip Battery pack(8800mAh)

(Option)

⑥ Charging Adapter(Belt clip battery pack type)

⑦ Cable(Belt clip battery pack)

⑧ 1.5x Loupe(Option)

⑨ 2.0x Loupe(Option)