



LEADING BRAND OF
**ELECTROSURGICAL
DEVICE**

EXPERIENCE CONVENIENT
ONE-STOP ELECTROSURGICAL
DEVICES PURCHASE



BAISHENG MEDICAL CO.,LTD.

CE0123 FDA 510(K) ISO 9001
13485


Surgery Solution Provider




Neurosurgery and
Brain Surgery



ENT and
Plastic Surgery



Vascular, Thyroid
and Breast Surgery




Cardiothoracic
Surgery




Hepatobiliary and
Pancreatic Surgery




Gastroenterology



Gynecology



Arthroscopy





OBS-100A

Monopolar Electrosurgical Unit
Max Power:100W

- 1.Portable and affordable
- 2.Suitable for general surgery and veterinary



OBS-100B

Cold Bipolar Coagulator
Max Power:80W

- 1.Low voltage output
- 2.Sparkless in work
- 3.Superb non-stick effect
- 4.Suitable for neurosurgery, ophthalmology, gynecology, and plastic surgery



OBS-100C

Dual cold RF Electrosurgical Unit
Max Power:120W

- Unique dual radio frequency design. Combine 4Mhz monopolar and 1.7Mhz bipolar output into one unit. Suitable for Plastic, Dermatology, ENT, Ophthalmology, Dentistry, Podiatry and Spinal surgery etc..



OBS-100C(I)

Intelligent Electrosurgical Unit
Max Power:100W

- 1.Portable in size
- 2.Available in both monopolar and bipolar modes
- 3.Convenient handle power control and memory function
- 4.Suitable for micro surgery, gynecology and arthroscopy



OBS-100C(II)

Cold RF Unit
Max Power:140W

- 1.Radio frequency
- 2.Precise in output and low temperature in work
- 3.No thermal tissue damage
- 4.Minimal collateral tissue damage
- 5.Suitable for gynecology, dermatology, ENT, ophthalmology, dentistry and podiatry



OBS-50

High Frequency Coagulator
Max Power:30W

- 1.Convenient handle power control
- 2.Precise coagulation
- 3.Suitable for dermatology, plastic surgery, proctology, gynecology, ophthalmology, urology, ENT, and dentistry



OBS-350A

Intelligent Versatile Electrosurgical Generator(Vessel Sealing)
Max Power:350W

- 1.Excellent compatibility
- 2.Available in vessel sealing
- 3.Suitable for general surgery, vascular surgery, hepatobiliary and pancreatic surgery, cardiothoracic surgery, gastroenterology, TURP and TUVF
- 4.Compatible with argon coagulation, cystoscopy, arthroscopy, laparoscopy and hysteroscopy



SS-100

Smoke Evacuation System
Max Power: 100W

Smoke evacuation system offers the effective way to smoke removal and protect operating room staff and patient from hazard.



OBS-350B SS-350

LCD Intelligent Electrosurgical Unit
Max Power:350W


High control precision of Argon. Intuitive graphic user interface, multifunctional, simple and practical in use. Vessel sealing function. Able to save users' custom parameters. Specialized power-feedback control technology.

Electrosurgical Pencil

OBS-Db(Disposable button control)

Applicable in general surgery


REF# 0039D stainless steel tip,sterilization pouch,1 pc/pouch,100 pouches/carton
REF# 0035D non-stick tip,sterilization pouch,1 pc/pouch,100 pouches/carton



OBS-Dr(Disposable rocker control)

Applicable in general surgery


REF# 0038D stainless steel tip,sterilization pouch,1 pc/pouch,100 pouches/carton
REF# 0030D non-stick tip,sterilization pouch,1 pc/pouch,100 pouches/carton




OBS-Ds(Disposable button control) Suction

Equip ed with smoke and liquid suction function.Flexible electrode tip can be extended from Suction Tube when cut or coagulate and shortened from Suction Tube when suck.Applicable in gastroenterology,hepatobiliary and pancreatic surgery.

REF# 0069S(-01.02.03) stainless steel tip,300mm length of medical tubing,blister,1 pc/blister,100 blisters/carton
REF# 0065S(-01.02.03) non-stick tip,300mm length of medical tubing,blister,1 pc/blister,100 blisters/carton
-01 equipped with 25mm length of suction tube,-02 with 85mm length of suction tube,-03 with 130mm length of suction tube.



Smoke evacuation



Applicable in gynecology & general surgery

REF# 0075S stainless steel tip,300mm length of medical tubing,blister,1 pc/blister,100 blisters/carton
REF# 0073S non-stick tip,300mm length of medical tubing,blister,1 pc/blister,100 blisters/carton



Smoke evacuation

Patent 200920313421.4

Applicable in gastroenterology, hepatobiliary and pancreatic surgery.

Flexible tip.Smoke elimination tube can move along with tip and equip ed with position lock.Replaceable with different tips



REF# 0079S stainless steel tip,300mm length of medical tubing, blister,1 pc/blister,100 blisters/carton

REF# 0071S non-stick tip,300mm length of medical tubing, blister,1 pc/blister,100 blisters/carton

OBS-Db(Disposable button control) Flexible

Applicable in gastroenterology and cardiothracic surgery. The tip can move to the ideal length,convenient in use.



REF# 0096S stainless steel tip,blister,1 pc/blister,100 blisters/carton
REF# 0095S non-stick tip,blister,1 pc/blister,100 blisters/carton

OBS-Df(Disposable foot control)

Applicable in general surgery

3.8mm



REF# 0048D stainless steel tip,sterilization pouch,1 pc/pouch,100 pouches/carton
REF# 0040D non-stick tip,sterilization pouch,1 pc/pouch,100 pouches/carton



OBS-Rb(Reusable button control)

Reusable with high temperature resistance,high pressure resistance and steam. Applicable in general surgery



REF# 0049D stainless steel tip,sterilization pouch,1 pc/pouch,100 pouches/carton
REF# 0045D non-stick tip,sterilization pouch,1 pc/pouch,100 pouches/carton

OBS-Da(Arthroscopic Ablator)


Incorporate a one-piece molded handle design and can be compatible with prevailing brands of electro-surgical generator.

Intended for ablation, tissue modification and coagulation of soft tissue in arthroscopic procedures such as hip, shoulder, ankle, wrist, elbow and knee

Features:

- 150mm length shaft in diameter 3.9mm
- Available in 45° and 90° angle head with operating area 4.9mm².
- Irrigation/suction device is available with a cannula and an adjustable flow control valve.
- Integrated, ergonomic handle design with hand controlled button operation
- Coagulation:50W(Coag)
- Ablation:200W(Cut)

REF No.Product Description

	#00201-45° ,45° angle flat brush tip #00201S-45° ,45° angle flat brush tip, with Suction		#00201A-90° ,90° angle round brush tip #00201AS-90° ,90° angle round brush tip, with Suction
	#00201-90° ,90° angle flat brush tip #00201S-90° ,90° angle flat brush tip, with Suction		#00202S-90° ,90° angle bloom tip, with Suction
	#00201A-45° , 45° angle round brush tip #00201AS-45° ,45° angle round brush tip, with Suction		#00202AS-90° ,90° angle lotus tip, with Suction

Package information: 25pcs per box,4 boxes per carton,100pcs per carton

Suction Coagulators

- 1.Suction and coagulation can be performed independently or simultaneously.
- 2.Increased efficiency during procedures requiring controlled fluid suction
- 3.Vented handle for suction control
- 4.Malleable cannula allows surgeon to match instruments to specific surgical need
- 5.Provide suction corrugated tubing and cleaning stylet

Control type:Button control
How to perform suction:place finger over suction vent on handle
Cannula diameter and length:4.2mm dia, 150mm L
Suction corrugated tubing length:300mm
Connecting cable and plug:3M with 3 pin plug
Not larger than standard coagulation power 30 watts is recommended for ENT operation.

REF	description	package details
0020S	Malleable cannula suction coagulator	Blister packing, 1 pc/blister, 25 bilsters/case, 4case/carton




Bipolar Irrigated Coagulator

Innovative cool HF electro-surgical device
REF# 0110



Monopolar Irrigated Knife

Innovative cool HF electro-surgical device
REF# 0100



Electrosurgical Electrode


Disposable LOOP tungsten electrode



- Disposable#0341 Loop shape 10X8 130mm
- Disposable#0342 Loop shape 10X10 130mm
- Disposable#0343 Loop shape 15X15 130mm
- Disposable#0344 Loop shape 20X20 130mm
- Disposable#0353 Loop shape 20X12 130mm
- Disposable#0354 Loop shape 20X15 130mm



- Disposable#0345 Square shape 10X5 130mm
- Disposable#0346 Square shape 10X10 130mm
- Disposable#0347 Square shape 15X10 130mm
- Disposable#0348 Square shape 20X10 130mm



- Disposable#0349 Triangle flag shape 15X10 130mm
- Disposable#0350 Triangle flag shape 17X12 130mm
- Disposable#0351 Triangle flag shape 20X15 130mm




- Disposable#0352 B shape 15X10 130mm


LOOPCUT II SERIES




- Disposable#0366 Square shape 10X5 79mm
- Disposable#0367 Square shape 15X10 80mm
- Disposable#0368 Square shape 20X10 82mm



- Disposable#0360 Loop shape 10X5 77mm
- Disposable#0361 Loop shape 10X10 78mm
- Disposable#0362 Loop shape 15X15 86mm
- Disposable#0363 Loop shape 20X12 85mm
- Disposable#0364 Loop shape 20X15 88mm
- Disposable#0365 Loop shape 20X20 92mm



- Disposable#0369 Triangle flag shape 15X10 82mm
- Disposable#0370 Triangle flag shape 20X15 87mm
- Disposable#0371 Triangle flag shape 17X12 85mm



- Disposable#0372 B shape 15X10 83mm

RELYPOINT SERIES



- Disposable #0306 Spatula Electrode 70mm
- #0306C Non-stick Spatula Electrode 70mm
- Reusable #0306(R) Spatula Electrode 70mm
- #0306C(R) Non-stick Spatula Electrode 70mm




Disposable #0312 Blade electrode 70mm
 #0312C Non-stick Blade electrode 70mm
 Reusable #0312(R) Blade electrode 70mm
 #0312C(R) Non-stick Blade electrode 70mm




Disposable #0324 Ø3mm Ball electrode 88mm
 #0324C Ø3mm Non-stick ball electrode 88mm




Disposable #0326 Blade electrode 100mm
 #0326C Non-stick blade electrode 100mm




Disposable #0327 Extended needle electrode 120mm
 #0327C Non-stick extended needle electrode 120mm



Disposable #0323 Ø5mm Ball electrode 89mm
 #0323C Ø5mm Non-stick ball electrode 89mm



Disposable #0325 Extended Spatula Electrode 89mm
 #0325C Extended Non-stick Spatula Electrode 89mm



Disposable #0308 Bent & Sharp Electrode 70mm



Disposable #0309 Bent & Blunt Electrode 70mm

MICROPOINT SERIES



Disposable #0317 Ø4mm Ball electrode 40mm #0317C Ø4mm Non-stick ball electrode 40mm
 Reusable #0317(R) Ø4mm Ball electrode 40mm #0317C(R) Ø4mm Non-stick ball electrode 40mm



Disposable #0319 Needle electrode 69mm #0319C Non-stick needle electrode 69mm
 Reusable #0319(R) Needle electrode 69mm #0319C(R) Non-stick needle electrode 69mm



Disposable #0310 Curved needle electrode 60mm #0310C Non-stick Curved needle electrode 60mm
 Reusable #0310(R) Curved needle electrode 60mm #0310C(R) Non-stick Curved needle electrode 60mm

Electrosurgical Electrode

MICROPOINT SERIES



Disposable #0315 Curved needle electrode 68mm #0315C Non-stick Curved needle electrode 68mm
 Reusable #0315(R) Curved needle electrode 68mm #0315C(R) Non-stick Curved needle electrode 68mm



Disposable #0313 Needle electrode 70mm #0313C Non-stick needle electrode 70mm
 Reusable #0313(R) Needle electrode 70mm #0313C(R) Non-stick needle electrode 70mm



Disposable #0311 Curved needle electrode 50mm #0311C Non-stick Curved needle electrode 50mm
 Reusable #0311(R) Curved needle electrode 50mm #0311C(R) Non-stick Curved needle electrode 50mm



Disposable #0318 Needle electrode 57mm #0318C Non-stick needle electrode 57mm
 Reusable #0318(R) Needle electrode 57mm #0318C(R) Non-stick needle electrode 57mm



Disposable #0321 Ø3mm Curved Ball electrode 61mm #0321C Ø3mm Non-stick Curved ball electrode 61mm
 Reusable #0321(R) Ø3mm Curved Ball electrode 61mm #0321C(R) Ø3mm Non-stick Curved ball electrode 61mm




Disposable #0301 Ø3mm Ball electrode 70mm #0301C Ø3mm Non-stick ball electrode 70mm
 Reusable #0301(R) Ø3mm Ball electrode 70mm #0301C(R) Ø3mm Non-stick ball electrode 70mm




Disposable #0302 Ø5mm Ball electrode 70mm #0302C Ø5mm Non-stick ball electrode 70mm
 Reusable #0302(R) Ø5mm Ball electrode 70mm #0302C(R) Ø5mm Non-stick ball electrode 70mm


ULTRALENGTH SERIES




Disposable #0303 Ø3mm Ball electrode 150mm
 #0303C Ø3mm Non-stick ball electrode 150mm




Disposable #0305 Ø5mm Ball electrode 150mm
 #0305C Ø5mm Non-stick ball electrode 150mm



Disposable #0316 Extended needle electrode 150mm
 #0316C Non-stick extended needle electrode 150mm








Disposable #0314 Extended Blade electrode 150mm
 #0314C Non-stick Extended blade electrode 150mm








Disposable #0307 Extended Spatula Electrode 150mm
 #0307C Extended Non-stick Spatula Electrode 150mm

RF SERIES



Electrosurgical Electrode







Full cover




Hydrogel over edge




No.#1030ac	Foam backing	Adult	Solid	Upright, With cable
No.#1031ac	Foam backing	Adult	Split	Upright, With cable
1 / , 25 / , 4 / .				

No.#1030a	Foam backing	Adult	Solid	Upright
No.#1031a	Foam backing	Adult	Split	Upright
1 / , 50 / , 8 / .				


Electrosurgical Pad




No.#1030b	Foam backing	Pediatric	Solid	Upright
No.#1031b	Foam backing	Pediatric	Split	Upright
1 / , 50 / , 8 / .				



No.#1030n	Foam backing	Neonatal	Solid	Upright
No.#1031n	Foam backing	Neonatal	Split	Upright
1 / , 50 / , 8 / .				




No.#1040a	Foam backing	Adult	Solid	horizon
No.#1041a	Foam backing	Adult	Split	horizon
1 / , 50 / , 8 / .				









No.#1040ac	Foam backing	Adult	Solid	horizon
No.#1041ac	Foam backing	Adult	Split	horizon
1 / , 25 / , 4 / .				


Safety design



①	No.#2030	Foam backing	Adult/Pediatric	Solid
	No.#2130	Non-woven fabric backing	Adult/Pediatric	Solid
②	No.#2031	Foam backing	Adult/Pediatric	Split
	No.#2131	Non-woven fabric backing	Adult/Pediatric	Split
1 / , 50 / , 4 / .				



Valleylab generator/Conmed generator

①	No.#4030	Foam backing	Adult	Split	Fit for Valleylab generator
	No.#4030C	Foam backing	Adult	Split, with cable	Fit for Valleylab generator
②	No.#4040	Non-woven fabric backing	Adult	Split	Fit for Conmed generator
	No.#4040C	Non-woven fabric backing	Adult	Split, with cable	Fit for Conmed generator
1  /  , 25  /  , 4  /  .					
















Reusable Silicon Electrosurgical Pad



①	#R1030a	Adult	Solid	1  / 
②	#R1031a	Adult	Split	
③	#R1030AC0801	Adult	Solid	
④	#R1030AC0806	Adult	Solid	

Connecting cables for patient plates

#R6060 (connect R1030a/R1031a)		3 meter length, For Alsa, Bard, Berchtold, Bove-Aaron 950/1200, Eschmann TD400/TD800, Martin, Meditom				
#R6070 (connect R1030a/R1031a)		3 meter length, For Valleylab				
#6060 (#0801-#0811)		3 meter length, For Alsa, Bard, Berchtold, Bove-Aaron 950/1200, Eschmann TD400/TD800, Martin, Meditom				
#6070 (#0803-#0811)		3 meter length, For Valleylab				
#6080 (#0804-#0811)		3 meter length, For Bowa, Conmed, Valleylab				
#6090 (#0806-#0811)		3 meter length, For Aspen, Bovie-Aaron 1250/2250/3250, Bowa, Conmed, Petas 500/600, Valleylab				
 #0801 Φ6.3mm plug	 #0802	 #0803	 #0804	 #0805	 #0806	 #0811

Electrosurgical Bipolar Forceps

OBS-Rp OBS-Dp

Overall length:100mm Tip size:0.5mm		Overall length:120mm Tip size:0.5mm		Overall length:110mm Tip size:0.5mm	
EU standard	Reusable #0550-01	EU standard	Reusable #0550-02	EU standard	Reusable #0550-03
US standard	Reusable #0560-01	US standard	Reusable #0560-02	US standard	Reusable #0560-03

Overall length:120mm Tip size:0.5mm		Overall length:110mm Tip size:0.5mm		Overall length:115mm Tip size:0.5mm	
EU standard	Reusable #0550-04	EU standard	Reusable #0550-05	EU standard	Reusable #0550-06
US standard	Reusable #0560-04	US standard	Reusable #0560-05	US standard	Reusable #0560-06

Overall length:170mm Tip size:0.5mm other tip size can be customized	EU standard	Reusable #0550-07
	US standard	Reusable #0560-07
Overall length:195mm Tip size:0.5、1.0mm other tip size can be customized	EU standard	Non-stick Reusable #0550-08 #0550-08(N)
	US standard	Non-stick Reusable #0560-08 #0560-08(N)
Overall length:160mm Tip size:0.5mm other tip size can be customized	EU standard	Reusable #0551-01
	US standard	Reusable #0561-01
Overall length:160mm Tip size:0.5mm other tip size can be customized	EU standard	Reusable #0551-02
	US standard	Reusable #0561-02
Overall length:180mm Tip size:0.5、1.0mm other tip size can be customized	EU standard	Reusable #0551-03
	US standard	Reusable #0561-03
Overall length:190mm Tip size:0.5、1.0mm other tip size can be customized	EU standard	Reusable #0551-04
	US standard	Reusable #0561-04

OBS-Rp OBS-Dp




Overall length:190mm Tip size:0.5mm other tip size can be customized 	EU standard		Reusable #0551-05
	US standard		Reusable #0561-05
Overall length:200mm Tip size:0.5, 1.0mm other tip size can be customized 	EU standard		Reusable #0551-06 Non-stick #0551-06(N)
	US standard		Reusable #0561-06 Non-stick #0561-06(N)

OBS-Ri OBS-Rin OBS-Di OBS-Din




Irrigation Forceps	Overall length:200mm Tip size:1.0mm		EU standard Reusable #0551-07
Irrigation Forceps	Overall length:195mm Tip size:1.0mm		EU standard Reusable #0551-08
Irrigation Forceps	Overall length:195mm Tip size:1.0mm		US standard Reusable #0561-07
Irrigation Forceps	Overall length:135mm Tip size:0.5mm		US standard Reusable #0560-09

OBS-Dp OBS-Dn

Disposable Bipolar Forceps Assembled With Cable

#0670	Overall length(exclude cable):195mm Tip size:1.0mm		
#0671	Overall length(exclude cable):200mm Tip size:1.0mm		

Disposable Irrigation Forceps

#0670i	Overall length(exclude cable):200mm Tip size:1.0mm		
#0671i	Overall length(exclude cable):225mm Tip size:1.0mm		

Bipolar Artery Forceps

Bipolar Artery Forceps No.#0542

#0542-01 Length 247mm
#0542-03 Length 193mm

#0542-02 Length 217mm
#0542-04 Length 166mm




#0542n-01 non-stick, Length 247mm
#0542n-03 non-stick, Length 193mm


#0542n-02 non-stick, Length 217mm
#0542n-04 non-stick, Length 166mm



Bipolar Artery Forceps No.#0543(Length 260mm)



Bipolar Artery Forceps No.#0543n(non-stick, Length 260mm)



Reusable Bipolar Forceps Cable

#0580 US Fitting		
#0581 EU Fitting		
#0582 EU Fitting		
#0584-1 EU Fitting		

Disposable bipolar forceps cable

#0580D	
--------	--

Monopolar laparoscopic electrode cable

REF # 08180810	
REF # 08180815	
REF # 08180817	
REF # 08180819	


OBS Disposable Skin Stapler

Model type:OBS-MW

Stapler Size/W-Type:7.0X4.0mm,Pin diameter 0.6mm

Staple No.:35pcs,25pcs,15pcs

Available specifications are 35W,25W,15W



7.0mm
4.0mm
W:7.0mmx4.0mm



Staple Remover

REF:B1

Removes atraumatically all brands of surgical skin staples;
Easy to handle making removal simple and safe.

Physical Therapy Electrode

Suitable for the use of TENS/EMS unit and other low or medium frequency stimulator.

Applied in neck,waist,legs and other parts of the body for the treatment of cervical syndrome,lumbar syndrome,bodybuilding and other physiotherapy

Available in two main series:electrode with stud and with lead wire


Highly conductive flim ensures low impedance and full dispersion of electrical current

Large selection of shapes and sizes maximizes the adaptability to multiple application location

Non-woven fabric backing provides ultimate flexibility and enhances patient comfort

High quality hydro gel proves to be totally biocompatible and not irritating

Durable material ensures its multiple-time uses without losing self-adhesive quality.




50X35mm



80X40mm



40X40mm



Φ 70mm



50X50mm



90X60mm



100X70mm

#9050 Electrode with stud

#9060



50X50mm



60X50mm

#9070



30X30mm

Electrode with cable



CE 0123 FDA 510(K)

ECG Monitoring Electrode (Suitable for Ordinary Monitoring, Anaesthesia & ICU)

301A




Φ 55mm Round Non-woven Fabric Backing Adult

301B



Φ 45mm Round Non-woven Fabric Backing Adult

301I radiolucenuy



Φ 45mm Round Non-woven Fabric Backing Adult

302B



Φ 45mm Peach Foam Backing Adult

302E



Φ 30mm Peach Foam Backing Child

302P



Φ 50mm Round Foam Backing Adult

Suitable for Stress Test & Holter

301F 35X45mm
Round Non-woven Fabric Backing Off-centered Adult

302F 35X45mm
Round Foam Backing Off-centered Adult

Patent No.:ZL200520063615.5




Disposable Subdermal Needle Electrodes


Single:
401-S-1000/13/0.4
402-S-1500/13/0.4
403-S-2000/13/0.4
404-S-2500/13/0.4

Twisted-pair:
405-W-1000/13/0.4
406-W-1500/13/0.4
407-W-2000/13/0.4
408-W-2500/13/0.4

Parallel-pair:
409-P-1000/13/0.4
410-P-1500/13/0.4
411-P-2000/13/0.4
412-P-2500/13/0.4




EU standard **US standard**



301G

Φ 30mm With Cable Non-woven Fabric Backing Newbom

Static Cardiogram Tab Electrode



301E

301E-01 22X34mm
301E-02 26X34mm
301E-03 28X34mm

Defibrillation Pad



OBS Defibrillation Pads are produced by high quality of hydrogel with excellent adhesivity through our unique automatical production line. It's hydrogel is overlapping design which can avoid the exposure of conductive film and provide homogeneous and uniform distribution of current on the surface to avoid burn happen.

OBS Defibrillation Pads are equipped with different connectors compatible with different brand of defibrillators.



▶ Defibrillation Pad Model: **OBS-DE/P**

The plug Model: **303A5030**

Connector compatible with: Philips Heartstart: FR2, FR2 (+), FR3, Heartstream Forerunner, Heartstart MRx, Heartstart XL, Heartstart XL+



▶ Defibrillation Pad Model: **OBS-DE/M**

The plug Model: **303A5035**

Connector compatible with:

Medtronic Physiocontrol, Osatu Bexen, Cardioline, Mindray



▶ Defibrillation Pad Model: **OBS-DE/W**

The plug Model: **303A3025**

Connector compatible with:

Nihon Kohden:
With pad adaptor JC-755V: TEC-5500 Series, TEC-7600 Series, TEC-7700 Series
With pad adaptor JC-855V: TEC-8300 Series



▶ Defibrillation Pad Model:

OBS-DE/Z1

The plug Model: **303A2545**

Connector compatible with:

Zoll Medical Corp



▶ Defibrillation Pad Model: **OBS-DE/Z2**

The plug Model: **303A6030**

Connector compatible with:

Zoll: E Series, M Series, R Series, AED PLUS, AED PRO



▶ The plug Model: **303A1525**

Connector compatible with:


CU Medical Systems, inc , Cmos Drake, Progetti



▶ The plug Model: **303A2060**

Connector compatible with:

Metrax Primedic





BAISHENG MEDICAL CO.,LTD.

Add: No.11 Fusheng Road, Xinhui District, Jiangmen, Guangdong, China

Tel.: +86-750-6628113 +86-750- 6691112

Fax.: +86-750-6616122

P.C.:529100

E_mail:baisheng@bs0750.com

Home page:http://www.bs0750.com

(Inter data)